F	NOTICE DE BELOTERO® VOLUME
GB	INSTRUCTIONS FOR USE FOR BELOTERO® VOLUME
D	GEBRAUCHSINFORMATION BELOTERO® VOLUME
I	ISTRUZIONI PER L'USO DI BELOTERO® VOLUME

F NOTICE DE BELOTERO® VOLUME

Description

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GB INSTRUCTIONS FOR USE FOR BELOTERO® VOLUME

Description

BELOTERO Volume is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

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BELOTERO Volume is presented in a single use pre-filled glass syringe, sterilized by moist heat. Each box contains one instruction leaflet, one or several syringe(s), two traceability labels per syringe and sterile CE-marked needle(s) and/or cannula(e) for single use only. The number of syringes, the dimensions and the number of needles and/or cannulae are stated on the external box.

Composition

Cross-linked sodium hyaluronate: 26.0 mg/ml

Phosphate buffer pH 7 q.s.: gel volume

The volume of the gel in each syringe is stated on the external box.

Indications

BELOTERO Volume is an injectable resorbable implant indicated to restore facial volumes, as for example to enhance the cheeks or the chin.

Posology and administration method

BELOTERO Volume is designed to be injected into the deep subcutaneous layers above the periosteum by a legally approved practitioner.

For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for restoring facial volume. A good knowledge of the anatomy and physiology of the site to be treated is required.

The treatment must be carried out under appropriate aseptic conditions.

BELOTERO Volume must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.

It is recommended to use one of the supplied needles.

To ensure optimal use of BELOTERO Volume it is recommended to assemble the needle according to the diagrams below.

If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle or the cannula.

BELOTERO Volume should be injected slowly. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. It is recommended not to inject more than 3 ml per treatment site at each session. Do not over-correct.

Gently massage the treated area after the injection to distribute the product uniformly.

Contra-indications

BELOTERO Volume is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate,
- In pregnant and breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.

Do not inject BELOTERO Volume into blood vessels.

Do not inject BELOTERO Volume into the labellar region.

Do not inject BELOTERO Volume into the tear-trough region.

Do not inject BELOTERO Volume into the lips.

Do not inject BELOTERO Volume into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).

Do not inject BELOTERO Volume for correction of superficial wrinkles and fine lines (injection of the product into the superficial dermis).

Do not inject BELOTERO Volume into an area previously treated with a permanent dermal filler.

Precautions for use

Before treatment the patient must be informed about the device, its contra-indications and possible side effects.

In the absence of available clinical data on tolerance and efficacy of the injection of BELOTERO Volume in patients with antecedents or with an active auto-immune disease or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Volume on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

BELOTERO Volume can be used to correct signs of facial fat loss (lipatrophy) in patients with human immunodeficiency virus.

It is recommended not to inject BELOTERO Volume in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or cheiloids.

It is recommended not to inject BELOTERO Volume into the periorbital region. BELOTERO Volume must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before completing the healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Volume must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Volume with the above-mentioned treatments.

No clinical data is available on the injection of BELOTERO Volume into an area already treated with another filling product.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Patients using anti-thrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of haematomas, nodules or bleeding at the injection site.

Do not transfer BELOTERO Volume into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe.

BELOTERO Volume must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.

Discard the syringe and the remaining product after use.

Do not re-sterilize and do not reuse due to the associated risks including infection.

The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

Incompatibilities

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Volume does not come into contact with such substances.

There is no known interaction with other local or loco-regional anaesthetics.

Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain sometimes accompanied by itching in the treated area. These reactions may last for a week.
- Haematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

• Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Assembly of needle to syringe

For optimal use of BELOTERO Volume it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.

2. Grasp the protective cap with the other hand and unscrew it.

3. Push & Twist the needle on the syringe until a resistance is felt. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.

4. Keep holding the Luer-lock and remove the sheath from the needle.

Storage

Store between 2 °C and 30 °C. Protect from light and freezing. Avoid mechanical shocks.

References

Updated documentation may be available from ANTEIS SA, Switzerland.

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Sterimedix Ltd